



G-4 Four-Mirror Glass High Mag

Primary Application:

Standard Gonio Lens for Static and Dynamic (no flange version) Gonioscopy

Features & Benefits:

- Crispest, clearest imaging with high magnification to appreciate details previously unattainable with Gonio lenses.
- All glass design provides superior clarity and durability compared to acrylic lenses.
- Four precisely angled mirrors provide full 360° views of the anterior chamber.
- Available in Flanged or No Flange (nf) contact options.

- No Flange (nf) contact design does not require interface solution allowing quick eye contact, reducing patient discomfort and facilitating the examination.
- No Flange version not rated for use with ophthalmic laser treatments.
- Available in large or small ring design to suit dexterity & personal preference.