

## HEINE BETA 200® Ophthalmoscope

### The ophthalmoscope with superior aspherical optics

The ophthalmoscope with unique aspherical optics and all the features of the BETA 200 S, except for the macula aperture/pinhole and the following lens range:

+in 1 D steps	1-10	15	20	40	D
-in 1 D steps	1-10	15	20	25	



- > Eliminates corneal reflex
- > Ideal for small pupils
- > Dustproof

<b>BETA 200 Ophthalmoscope</b>	2.5 V	3.5 V
BETA 200 Ophthalmoscope	C-01.30.100	C-02.30.100
BETA 200 M2 Ophthalmoscope, like BETA 200, except with fixation star aperture without polar co-ordinates and with pinhole instead of small spot	C-01.30.102	C-02.30.102
XHL Xenon Halogen spare bulb	X-01.88.069	X-02.88.070

#### BETA 200 apertures



Slit, fixation star with polar co-ordinates, cobalt blue filter, large spot, small spot, hemispherical spot



With red-free filter

#### BETA 200 M2 apertures



Slit, fixation star, cobalt blue filter, large spot, pinhole, hemispherical spot



With red-free filter