

## Bish



### op Alignment

### Prism Set

(11311)

The Bishop Alignment Set contains a metal plate, magnetic prisms of the following sizes: 0.5, 1, 2, 4, 6, 8, 10, 12, 15, 20, 25, 30, 35, 40, 45, 50 Diopter, and a red lens. The magnetic prisms in the Bishop set have magnets on both the x and y axis of the prism. The Bishop tool was developed to provide advantages in **measuring large angle strabismus** prior to correction. The Bishop Alignment Prism Set also comes with a galvanized bishop plate that allows multiple prisms to be fixated upon it. Click the link (left) to learn more.